

# WizPure™ PCR FDMix

(Aliquot & lyophilized in 8-strip PCR tube)

- W1601 96 rxn
- W1601-5 480 rxn

## Description

WizPure™ PCR FDMix combines all the reagents necessary for successful routine PCR in a convenient individually aliquot and lyophilized in single tube. WizPure™ PCR FDMix is an economical, highly efficient and ready-to-use PCR premix which can amplify templates of up to 5 kb. The PCR FDMix contains Taq DNA Polymerase, dNTPs, MgCl<sub>2</sub>, enhancer and stabilizer.

The amplification products are compatible with TA cloning. WizPure™ PCR FDMix yields excellent and consistent results in routine PCR reactions as well as high-throughput PCR genotyping, colony PCR, RT-PCR and PCR cloning. PCR reactions can be directly loaded onto an agarose gel without the additional need of loading buffer and dyes.

WizPure™ PCR FDMix without gel loading dyes is also available for applications when loading dyes are undesired.

## Kit Contents

Contents	W1601	W1601-5
WizPure™ PCR FDMix	96 tubes	480 tubes

## Applications

WizPure™ PCR FDMix is suitable and tested for amplification of genomic targets ranging from 100 bp to 4 kb and of episomal targets (lambda phage; plasmids) up to 10 kb under various reaction conditions.

- High through-put PCR
- Routine diagnostic PCR requiring high reproducibility
- DNA sequencing template preparation

## Storage Conditions

Upon receipt, store all components at -20°C.

## Note

Do not contaminate the WizPure™ PCR FDMix with primers and template DNA used in individual reactions. Thaw and mix all components thoroughly, spin down shortly and chill on ice.

## Quality Control

No endonuclease activity, nicking activity, exonuclease activity, or priming activity has been detected.

Quality Authorized by : Jamie Ahn



## Protocol

This standard protocol applies to a single reaction where only template, primers, and water need to be added to the PCR FDMix. For multiple reactions, scale-up volume of reaction components proportionally. All reagents should be thawed on ice, gently mixed and briefly centrifuged before use.

1. Prepare the reaction mixture as following table.

### Reaction Conditions

Component	20 µl reaction	Final Conc.
PCR FDMix	1 tube	1X
10µM Forward Primer	0.2 ~ 2.0 µl	0.1~1.0 µM
10µM Reverse Primer	0.2 ~ 2.0 µl	0.1~1.0 µM
Template DNA	1 ~ 5 µl	< 250 ng
Water, RNase-Free	up to 20 µl	

**NOTE:** In general, use greater than 0.5 µM primers for sensitivity and less than 0.5 µM for specificity.

**NOTE:** Recommended amount of template per PCR reaction:

- < 50 ng plasmid or
- < 500~1000 ng genomic DNA or
- 2 µl of a 100 µl single plaque eluate or
- one single bacterial colony

2. Ensure reactions are mixed thoroughly by pipetting or gentle vortexing followed by a brief spin in a microcentrifuge.  
(Optional) Overlay reactions with one-half volume PCR-grade mineral oil when not using heated lid on thermal cycler.
3. Transfer tubes into a PCR instrument and run as following table.

### PCR Conditions

Step	Temp (°C)	Time	Cycle
Initial Denaturation	95	5 min.	1
Denature	95	10 ~ 60 sec.	25 ~ 40
Anneal	50~65	10 ~ 60 sec.	
Extend	72	60 sec./kb	
Final Extension	72	5 min.	1

**NOTE:** Cycling conditions may need to be optimized, depending on different primer and template combinations. For example, raise the annealing temperature to prevent non-specific primer binding, increase extension time to generate longer PCR products.

4. After cycling, maintain the reactions at 4°C or store at -20°C until ready for analysis.

**RUO** Research Use Only